
INSTRUCTIONS FOR USE

Isiris™ α

- GB** STENT REMOVAL DEVICE
- FR** DISPOSITIF DE RETRAIT DE SONDE
- DE** VORRICHTUNG ZUR URETERSCHIENEN-ENTFERNUNG
- IT** DISPOSITIVO DI RIMOZIONE DELLO STENT
- ES** DISPOSITIVO DE RETIRADA DE TUTOR
- PT** DISPOSITIVO DE REMOÇÃO DE STENTS
- NL** STENTVERWIJDERINGSINSTRUMENT
- GR** ΣΥΣΚΕΥΗ ΑΦΑΙΡΕΣΗΣ STENT
- RU** УСТРОЙСТВО ДЛЯ УДАЛЕНИЯ СТЕНТА
- PL** URZĄDZENIE DO USUWANIA STENTU
- CZ** PROSTŘEDEK PRO ODSTRANĚNÍ STENTU
- TR** STENT ÇIKARMA CİHAZI
- FI** STENTIN POISTOLAITE
- SE** ENHET FÖR STENTBORTTAGNING
- DK** STENTFJERNELSESANORDNING
- NO** ENHET FOR FJERNING AV STENT
- CN** 支架移除器械
- BG** ИЗДЕЛИЕ ЗА ИЗВАЖДАНЕ НА СТЕНТ
- HU** SZTENTELTÁVOLÍTÓ ESZKÖZ
- SK** POMŮCKA NA VYBERANIE STENTU
- RO** DISPOZITIV DE ÎNDEPĂRTARE A STENTULUI
- RS** UREDAJ ZA UKLANJANJE STENTA
- 阿拉伯文** ةمراع دلا ةل ازا زاهج

Manufactured by:
Ambu A/S
Baltorpbakken 13
2750 Ballerup
Denmark

Distributed by:
Coloplast A/S
Holteham 1
3050 Humlebaek
Denmark

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Hongkong – Keymax Technology, Unit 2103, 21/F., Kodak House II, 39 Healthy Street East, North Point- Hong Kong.
Phone :+852 31508913 / Fax :+852 31508412 – info@keymax.com.hk.

Hungary – C.P.P Budapest Kft, Amerikai UT 33, H-1145 Budapest.
Phone : +36 14223700 // Fax: +36 1 422 37 08 – eva.kandiko@cpcbudapest.hu.

Iran – Bartar Co, Bldg N°33 Sarparast ST Cross, Bozorgmehr St, Sth Felestin, Tehran.
Phone : + 98 21664 87360 – import1@bartarmedical.com.

Iraq – Al Zahrawi Scientific Bureau, D.605 ST. 13 H.7/A, Al Mansur, Baghdad.
Phone : +964 770 271 6356 – mayad@zahrawiscientific.com.

Israel – Sela Medical, 14 Imbar St., Petach Tikva 49511.
Phone : +972 3 9043433 / Fax: +972 3 9045474 – maya@selamedical.co.il.

Jordan – Al Eiman Drug store, Al Bayader, Industrial St., Port Street Beside the Civil Defence, 11810 Amman.
Phone : +962 6 581 0108 / Fax: +962 6 581 8850 – h.saifi@eiman-eds.com.

Korea – Suntek Medical Systems , 3rd Fl. Zero BD- 84 Gil, 14 Teheran-Ro , Kangnam-Ku-Seoul, 06178.
Phone : +82 2 5143214 / Fax +82 2 547 2446 – cypark99@chol.com.

Kuwait – Mohamed Naser Al-Hajery & Sons Co Ltd, street 67, block 73, Shuwaikh Industrial Area, Shuwaikh.
Phone : +965 248 31000 / Fax: +965 248 34029 – ahmed.saleem@hajery.com.

Lebanon – Saramed SARL, Liberty Building, Kobayyat St., Rmeil,2079 1204 – Beirut.
Phone : +961 1577 111 / Fax: +961 1575 600 – m.hasswany@saramed.com.

Libya – Alemtiyaz Al Awal, Alfuwahat Dubai Street, Building number 20, Benghazi.
Phone : + 218 91 9999 654 – youssef@alemtiyazalawal.com.

Lithuania – Skirgesa Ltd, Ekskavatorininku str. 1b, LT-52461 Kaunas.
Phone : +370 37 45 77 46, 45 21 68 / +370 37 45 81 61 – chirurgija@skirgesa.lt.

Malaysia – Pro Element Technologies , 18-34-A3 Persiaran Gurney, Gurney Tower-10250 Penang.
Phone : +604 226 32 86 / Fax +604 227 3286 – chiew@pro-element.com.my.

Malta – Pharma-Cos Ltd., Pharma-Cos House Triq. C. Portanier, 1732 Santa Venera HMR 11.
Phone : +356 2144 1977 / Fax: +356 2144 1388 – mmifs@pharma-cos.com.

Mexico – Geoworks Durango 225-204, Col. Roma – Mexico, D.F 06700.
Phone : +52 55 5514 0381 – geoworks@geoworks.com.mx.

Morocco – Sahmed, 2, rue des Acacias & Boulevard Y El Mansour El Maarif, 20100 Casablanca.
Phone : +212 2298 1177 / Fax: +212 2298 1178 – logistics@claesmedical.com.

Norway – Endomed A.S., Prof. Birkenlands vei 34B – 1081 Oslo.
Phone : +47 23 17 24 00 / Fax: +47 23 17 24 01 – vibeke.birkeli@endomed.no.

Oman – Etihad Medical Services LLC, Post Box 176, Postal Code 103, Muscat.
Phone : +968 246 16085 – rahul.ramakrishnan@etihadmedical.com.

Pakistan – FW. Distributors, FW Building, Opposite Poonch House, Adamjee Road, Saddar – Rawalpindi.
Phone : +92 515516729 / Fax: +92 515120061 – fw@fwdpak.com.pk.

Paraguay – Albatros S.A.C.I., Dr. N. Benitez 404, Asuncion.
Phone : +595 21 441 407 – mthiede@albatros.com.py.

Peru – Famedical EIRL – Calle Emilio Garcia Rosell 773, Urb. Las Brisas – Lima.
Phone : +511 723 26 15 – acomercial@famedical.com.pe.

Pologne – Skamex, Spolka z ogr. odp. Sp. k., Czeszochowska 38/52, 93-121 Łódź.

Phone : +48 42 677 14 11 / Fax: +48 42 672 40 10 – Katarzyna.kowalczyk@skamex.com.pl.

Qatar – Al Danah Medical Company W.L.L, Gate N°4 Naser Bin Khaled Complex, Salwa Rd , PO Box 14485, Doha

Phone : +974 4469 1122 / Fax : +974 4469 1124 – info@aldanahmedical.com.

Romania – Steril Romania, 3-5, Metalurgiei Blvd., Sector 4 041831, Bucharest.

Phone : +4021 460 07 79 / Fax: +4021 460 00 14 – import@steril.ro.

Russia – MPA Medical Partners, 8th Martha Str, 1 bld.12- Moscow 127083.

Phone : +7-495 921 30 88 / Nkibko@mpamed.ru.

Russia – UROMED, Vvdenskogo str., 1 bld. 1 – Moscow 117342.

Phone : +7-495 783 68 11 / forsova@uromed-m.ru.

Saudi Arabia – JAT, 80 Amr Ibn Al-Aas street, Al Wisham Area Office 104, Riyadh.

Phone : +966 114040470 / Fax : +966 11 4010041 – maman@janatarab.com.

Serbia – Magna Pharmacia, Milutina Milankovica 7b, Novi Beograd – Belgrade 11070.

Phone : +381 11 36 22 500 / Fax : +381 11 36 22 600 – bojan.peric@magnapharmacia.rs.

Singapore – Mediconsult Asia Pte Ltd, 2 Ang Mo Kio St. 64, Nr 04-03C, Econ Building Singapore – 569084

Phone : +65 6483 2370 – molly@mediconsultasia.com.sg.

Slovakia – SEQ-MED S.R.O., Biskupicka 21 , 986 01 Filakova.

Phone : +421.2 555 766 05 / Fax : +421 47 438 25 60 – info@seqmed.sk.

Slovenia – Apollonia d.o.o., Zaloska 155, 1000 Ljubljana.

Phone : +386 1 548 35 10 / Fax: +386 1 540 03 97 – apollonia@amis.net.

South Africa – Nu Angle, Postnet Suite 93, Private Bag X33, Rivonia.

Phone : +27 114674145 – dawn@nuangle.co.za.

Taiwan – LEGENE ENTERPRISE CO. Ltd. 6f, No.13, Sec. 2nd., Zhonghua Rd. Xinzhuang Dist. New Taipei City 242.

Phone : +02 2997 0186 / Fax : +02 2997 0336 – db1727@yahoo.com.tw.

Thailand – Ouiheng International Healthcare CO Ltd, 46, 46/2 Charansanitwongs 40 Bangyikhan,

Bangplad – 10700 Bangkok.

Phone : +66 2-886-8600 – sansanee@ohi-hc.com.

Tunisia – Medico Chirurgical, 31 BD Hedi Saidi , 1005 Tunis.

Phone : +216 71899 750 / Fax: +216 7189 7050 – medico.import@gnet.tn.

Turkey – Çapa Medikal, Yassioiren Mh. Agackakan Sk. No : 15A, Arnavutkoy, 34277 Istanbul.

Phone : +90 2021 444 8717 / Fax : +90 2021 671 9070 – cagdas.boyrac@capamedikal.com.tr

United Arab-Emirates – Pharmatrade, P.O. Box 11397, Dubai.

Phone : +971 4567 6135 / Fax: +971 4440 5101 – mzainaldeen@pharmatradeuae.com.

Uruguay – H. Strattner –Neomed SRL ,Francisco Llambi 1556, Montevideo, Uruguay.

Phone : +598 2 622 8326/ Fax: +598 2 622 5567 – rocioabella@neomed.com.uy.


USA – Coloplast Corp. Minneapolis, MN 55411.

Phone : +1-612-337-7800 – www.coloplast.com. For customer service, please call 1-800-258-3476 in USA; Outside USA, contact your local Coloplast representative.

Vietnam – Toan Anh CO. LTD, 162 Co Giang St., Co Giang Ward, Dist.1, Hochiminh City.

Phone : +848 377 55569 – binh@toananhcompany.com.

Isiris α

The Coloplast logo is a registered trademark of Coloplast A/S. © [2018-01.]
All rights reserved. Coloplast A/S, 3050 Humlebaek, Denmark.

Content	Page
Instructions for use	7 – 19
Mode d'emploi	20 – 32
Bedienungsanleitung	35 – 45
Istruzioni per l'uso	46 – 58
Instrucciones de uso	59 – 71
Instruções de utilização	72 – 84
Gebbruksaanwijzing	85 – 97
Οδηγίες χρήσης	98 – 110
Инструкция по применению	111 – 123
Instrukcja obsługi	124 – 136
Návod k použití	137 – 149
Kullanım talimatları	150 – 162
Käyttöohjeet	163 – 175
Bruksanvisning	176 – 188
Brugsanvisning	189 – 201
Bruksanvisning	202 – 213
使用说明书	214 – 225
Инструкции за употреба	226 – 238
Használati útmutató	239 – 251
Návod na použitie	252 – 264
Instrucțiuni de utilizare	265 – 277
Uputstva za upotrebu	278 – 290
304 – 291	مراذخ ست الال تا مری لعت

ISIRIS™ α STENT REMOVAL DEVICE

Contents	Page
1. Important information – Read before use	8
2. System description	10
3. Technical product specifications	12
4. Inspection and preparation of Isiris α	14
5. Operating Isiris α	15
6. After use.....	17
7. Trouble shooting	17
Appendix 1: Electromagnetic compatibility.....	19
Appendix 2: Warranty	19

1. Important information – Read before use

NOTE

Read these safety instructions carefully before using the Isiris Stent Removal Device. The Instructions for use may be updated without further notice. Copies of the current version are available upon request.

1.1. Instructions

The Isiris System is for use in hospital or in a medical office environment.

Isiris α is a sterile, single-use device designed for patients, accepting a scope diameter of 5.5 mm, implanted with a double loop ureteral stent, ranging from size CH4.8 to CH9.

Please be aware that these instructions do not explain or discuss clinical procedures. They describe only the basic operations and precautions related to the operation of the Isiris α . Before initial use of the Isiris α , it is essential for operators to have received sufficient training in clinical cystoscopic techniques and to be familiar with the intended use, warnings, cautions, notes, indications and contraindications mentioned in these instructions.

1.2. Intended use

Isiris α is a sterile, single use flexible cystoscope designed for removal of double loop ureteral stents accessible in the bladder via an urethral insertion in adults. Isiris α has been designed to be used with the reusable Isiris Monitor to visualize the observations obtained by Isiris α .

1.3. Contraindications

- Febrile patients with urinary tract infections (UTIs) or severe coagulopathy.
- Patients with acute infection (acute urethritis, acute prostatitis, acute epididymitis).
- Patients with known unpassable urethral stricture.

1.4. Warnings, cautions and notes

Appropriate warnings, cautions and notes are given describing potential safety hazards associated with the use of the Isiris α . The information given in this instruction serves only to instruct in the correct handling of the device.

WARNING

1. Only to be used by trained physicians, trained in clinical endoscopic techniques and procedures as incorrect use can cause harm to the patient.
2. Images produced by the Isiris System must not be used as an independent diagnostic tool of any pathology. Health care professionals must interpret and substantiate any finding by other means and in the light of the patient's clinical characteristics to avoid wrong diagnosis.
3. Do not use the Isiris System for any other purpose than its intended use as it may lead to patient harm.
4. Discard the product if the Inspection and Preparation (section 4) of the product fails as it can cause infection or patient injury.
5. While withdrawing the scope, the distal tip must be in a neutral and non-deflected position. Do not operate the bending lever, as this may result in injury to the patient.
6. Do not attempt to clean and reuse the Isiris α as it is a single use device. Reuse of the product can cause contamination leading to infections, major allergic reactions or patient injury.
7. The Isiris System is neither MRI safe nor MRI compatible.

8. Do not use the Isiris System during defibrillation as this may result in electrical shock to the user.
9. To avoid risk of electric shock, this equipment must only be connected to a supply mains with protective earth. To disconnect the Isiris System from mains, remove the mains plug from the wall outlet.
10. Check that the grasper is not protruding from the endoscope irrigation channel prior to insertion. The grasper teeth are sharp and can cause damage.
11. Excessive force should never be used when operating Isiris α , as doing so may result in mucosal injury.
12. Always watch the live endoscopic image on the display and adequately monitor the patient when inserting or withdrawing the scope, and when operating the grasper, bending section or suctioning. Failure to do so may result in patient injury.
13. Electronic equipment may affect the normal function of the Isiris System. If the Isiris System is used adjacent to or stacked with other equipment, observe and verify normal operation of the Isiris System prior to using it. The Isiris System may cause interference or disrupt equipment operations nearby. It may be necessary to adapt procedures for mitigation, such as reorientation or relocation of the equipment.
14. Do not activate the grasper when the distal end is inside the urethra as this may result in serious injury to the patient.
15. Do not enter any part of the Isiris α into the ureter as this may result in injury to the patient.
16. The Isiris α is not designed for introduction of endoscopic instruments. Do not try to introduce any device or contrast agent through the irrigation connector, as this may damage the graspers functionality.
17. The maximum allowed irrigation pressure limit is 200 cmH₂O or 0,2bar. Exceeding this limit may result in serious injury to the patient.
18. The distal end of the endoscope may get warm due to heat from the light source. Avoid long periods of contact between the tip of the device and the mucosal membrane as long, sustained contact with the mucosal membrane may cause mucosal injury.
19. During removal of the double loop stent, it is the user's responsibility to evaluate the pull strength and adapt the removal force and speed. Never exceed the usual pulling force as it may damage the patient's urinary tract.
20. Always perform a visual check according to the instructions in section 4 before placing the Isiris α in a waste container to ensure that no parts are left inside the patient.
21. Isiris α is considered infected after use and must be disposed of in accordance with local guidelines for collection of infected medical devices with electronic components.
22. Disconnect Isiris Monitor from any mains power supply, remove any accessories and make sure Isiris Monitor is completely turned off before cleaning and disinfection. The equipment is disconnected from mains by removing the mains plug from the wall outlet

CAUTION

1. Do not use sharp instruments to open the pouch or cardboard box as this could damage the sterile barrier and compromise the scope. Use standard medical practice for handling sterile equipment.
2. Make sure to have a suitable backup system readily available in case a malfunction should occur.
3. Keep the Isiris α handle dry during preparation and use, as liquid can cause a short circuit and affect the displayed image.
4. The Isiris System consists of the parts described in section 2. Failure to comply with this may reduce safety and efficiency.

5. Be careful when handling the distal tip of the insertion cord and do not allow it to strike other objects, as this may result in damage to the equipment. The lens surface of the distal tip is fragile and visual distortion may occur.
6. The scope connector must not be twisted when inserted into the monitor plug - this could potentially damage the scope connector or the monitor plug.
7. Do not exert excessive force on the bending section, as this may result in damage to the equipment. Examples of inappropriate handling of the bending section include:
 - Manual manipulation and bending such as performing retrograde vision inside bladder.
 - Operating it in any case where resistance is felt.
8. Be sure to connect the irrigation line properly to the irrigation connector before introducing fluid. Failure to do so may result in the fluid spilling from the irrigation connector, which may harm the device.
9. Secure syringes properly on the irrigation connector before aspiration is applied to avoid leakage.
10. This device contains the following substance defined as CMR 1B in a concentration above 0.1 % weight by weight: Cobalt; (CAS No. 7440-48-4; EC No. 231-158-0). Current scientific evidence supports that medical devices manufactured from stainless steel alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects for all patient populations.
11. US federal law restricts this device to sale by or on the order of a licensed health care practitioner.

NOTES


- The manufacture is not responsible for any damage to the system or patient resulting from incorrect use.

2. System description


The Isiris α can be connected to the Isiris Monitor. For information about the Isiris Monitor, please refer to the Isiris Monitor Instructions for use.

2.1. System parts


Isiris α


The Isiris System consists of Isiris α and Isiris Monitor. Isiris α is a sterile single use device. Isiris Monitor is reusable.

Isiris a**Protective sheath**

Protects the insertion cord during transport and storage. Remove before use

**2.2. Explanation of symbols used**

Symbols for Isiris a	Indication
	Working length of Isiris a insertion cord
	Maximum insertion portion width (Maximum outer diameter)
	Maximum length of protruded open grasper and minimum grasper opening distance
	Field of view
	Connection for Isiris a
	Electrical Safety Type BF Applied Part
	Single sterile barrier system. Sterilized using ethylene oxide
	UL Recognized Component Mark for Canada and the United States

Symbols for Isiris a	Indication
 CAS No. 7440-48-4	Contains hazardous substances
	Global Trade Item Number
	Medical Device
	Country of manufacturer
	Russian mandatory [Gosstandard's mark of conformity]
	CAUTION, Federal (USA) Law restricts this device to sale by or on the order of a physician
	Temperature limit
	France and Germany : Symbol of contribution for sorting and recycling the packaging materials
	Do not use if package is damaged and consult instruction for use
	Relative humidity limitation
	Atmospheric pressure limitation
	Recycling
	Authorised representative in Switzerland

3. Technical product specifications

3.1. Isiris a specifications

Optical system

Field of View	85°
Depth of Field	8 – 19 mm
Illumination method	LED

Insertion portion

Deflecting section ³	90° down, 80° up
Insertion cord diameter	5.0 mm (15 Fr/0.20 inch)
Distal end diameter	5.4 mm (16 Fr/0.21 inch)
Maximum diameter of insertion portion	5.5 mm (16.5 Fr/0.22 inch)
Working length	390 mm (15.4 inches)

³ Please be aware that the deflection angle can be affected if the insertion cord is not kept straight

Grasper

Maximum length of protruded grasper	18 mm (0.71 inch)
Minimum grasper opening distance	4.5 mm (0.16 inch)

Irrigation connector

Connection type	Standard Luer connector
-----------------	-------------------------

Operating environment

Temperature	10 – 40 °C (50 – 104 °F)
Relative humidity	30 – 85 %
Atmospheric pressure	80 – 106 kPa
Altitude	≤ 2000 m

Storage

Temperature	10 – 25 °C (50 – 77 °F)
Relative humidity	10 – 85 %
Atmospheric pressure	50 – 106 kPa

Transportation

Temperature	-10 – 55 °C (14 – 131 °F)
Relative humidity	10 – 95 %
Atmospheric pressure	50 – 106 kPa

Sterilisation


Method of sterilisation	EO
-------------------------	----

4. Inspection and preparation of Isiris α

4.1. Visual inspection of Isiris α

1. Check that the pouch seal is intact.
2. Make sure to remove the protective sheath from the insertion cord.
3. Check that there are no impurities on the product.
4. Check that there is no evidence of shipping damage or other damage such as rough surfaces, sharp edges or protrusions which may harm the patient.

4.2. Inspection of the Image

1. Connect Isiris α to Isiris Monitor by plugging the white connector with red arrow into the corresponding red port on Isiris Monitor.
2. The connector must not be twisted when inserted into the plug - this could potentially damage the scope connector or the monitor plug.
3. Verify that a live video image appears on the screen.
4. Point the distal end of Isiris α towards an object, e.g. the palm of your hand.
5. Adjust the image brightness and contrast preferences on Isiris Monitor if necessary – please refer to the Isiris Monitor *Instructions for use* for details. It is not possible to adjust the white balance since this is predefined in the system.
6. If the object cannot be seen clearly, wipe the lens at the distal end using a sterile cloth.

4.3. Preparation of Isiris α

1. Carefully slide the deflection lever up and down in each direction until it stops. Confirm that the deflective tip functions smoothly and correctly.
2. Slide the deflection lever slowly to its neutral position. Confirm that the deflective tip returns smoothly to a neutral position.
3. Activate the grasper by pressing the button and observe that the opening, closing and grasper retraction movements are completely executed. Perform this check in a straight position and in a full up and down deflection.
4. Connect to the irrigation line, turn on irrigation and make sure that irrigation flushes through the scope.
5. Aspiration might be required during the procedure. Make sure to prepare a 50 cc syringe for this. When required, adapt the syringe to Isiris α and apply an aspiration force according to the wanted effect. Disconnect, empty the syringe and reconnect it should there be a larger quantity of liquid to aspirate.

4.4. Additional preparations

In addition to the Isiris α , prepare the following:


- Medical grade water based lubricant
- Irrigation line with Luer connector and irrigation bag
- 50 cc syringe

5. Operating Isiris a

5.1. Holding Isiris a


The hand piece of Isiris a can be held in either hand. Use the thumb to move the deflection lever and the index finger to operate grasper activation button. The hand that is not holding Isiris a can be used to advance the insertion cord into the patient's urinary tract.

5.2. Manipulating the tip of Isiris a


The deflection lever is used to flex and extend the tip of Isiris a up and down. Moving the deflection lever downward will make the tip deflect up. Moving it upward will make the tip deflect down (European movement).

The insertion cord should be held as straight as possible at all times in order to secure an optimal tip deflection angle and grasper functionality.

5.3. Irrigation with saline solution

Connect the irrigation line Luer connector to Isiris a irrigation connector and make sure irrigation flushes through the irrigation channel. Check that the Luer connector is correctly locked. Failure to do so may result in the saline solution spilling, which may harm the device and reduce the irrigation flow.

5.4. Aspiration

Aspiration might be required during the procedure. Make sure to prepare a 50 cc syringe for this.

When required, adapt the syringe to Isiris α and apply an aspiration force according to the wanted effect; disconnect, empty the syringe and reconnect it should there be a larger quantity of liquid to aspirate.

If needed, partly or completely empty the bladder to the wanted volume by using a syringe adapted to the irrigation connector


Apply a vacuum of 0.9 bar or less when aspirating.

5.5. Grasper

The grasper is integrated in Isiris α and cannot be removed or detached. It is activated by firmly pressing the grasper activation button as follows:

1. Grasper opening: press the grasper button approximately half-way down.
2. Grasper closing over the double loop stent: continue pressing the grasper button to the maximum possible.
3. Release of the double loop ureteral stent from the grasper: at any time the double loop ureteral stent can be released by releasing the button. In rare cases, a second pressing and release of the button may be necessary.

Catching the double loop ureteral stent is performed the following way and always under visual control: once the double loop ureteral stent is in the field of view of Isiris α , open the grasper as indicated above, advance Isiris α with the opened grasper until getting in contact with the double loop ureteral stent, close the grasper as indicated above and maintain the button fully pressed. Release any deflection of the deflective tip. The removal of the double loop ureteral stent can start by gently pulling Isiris α .

5.6. Insertion of Isiris α

Lubricate the insertion cord with a medical grade water based lubricant to ensure the lowest possible friction when Isiris α is inserted into the patient's urethra. It is usual practice to lubricate the patient's urethra using the same water based lubricant.

If the camera image of Isiris α becomes unclear the tip can be cleaned by removing Isiris α and cleaning the tip with a piece of sterile gauze. Continue the procedure until a satisfactory image is obtained.

5.7. Withdrawal of Isiris α

When withdrawing Isiris α , make sure that the deflection lever is in the neutral position. Slowly withdraw Isiris α while watching the live image on Isiris Monitor.

By the end of procedure, put the distal tip in its neutral and non-angled position and slowly withdraw the Isiris α without touching the bending lever. Always observe the live endoscopic image while operating the grasper and manoeuvring and withdrawing the Isiris α .

6. After use

Visual check

1. Are there any missing parts on the deflective tip, lens, insertion cord, or grasper? If yes, then take corrective action to locate the missing part.
2. Is there any evidence of damage on the deflective tip, lens, insertion cord, or grasper? If yes, then examine the integrity of the product and conclude if there are any missing parts.
3. Are there cuts, holes, sagging, swelling or other irregularities on the deflective tip, lens, insertion cord, or grasper? If yes, then examine the product to conclude if there are any missing parts. In case of corrective actions needed (step 1 to 3) act according to local hospital procedures. The elements of the insertion cord and grasper are radiopaque.

Final steps

Disconnect Isiris α from Isiris Monitor and dispose of Isiris α in accordance with local guidelines for collection of infected medical devices with electronic components. Switch off Isiris Monitor by pressing the power button for at least 2 seconds.

7. Trouble shooting

If problems occur with the Isiris α , please use this trouble-shooting guide to identify the cause and correct the error.

Absent or reduced aspiration capability.

Cause	Action
Irrigation channel blocked	Clean the irrigation channel by flushing the irrigation channel with sterile saline using a syringe.

No live image on the left side of the screen but User Interface is present on the display.

Cause	Action
Isiris a is not connected to Isiris Monitor.	Connect Isiris a to the red port on Isiris Monitor.
Isiris Monitor is operated in yellow file management tab or green settings tab.	Return to live image by pressing the blue Live image tab or the red button with a black X in the top right corner.
Isiris Monitor and Isiris a have communication problems.	Restart Isiris Monitor by pressing the ON/OFF button for at least 2 seconds. When Isiris Monitor is OFF restart by pressing ON/OFF button once more.
Isiris a is damaged.	Replace Isiris a with a new one.

The image shown to the left is frozen.

Cause	Action
A communication error has occurred between Isiris a and Isiris Monitor.	Restart the system by pressing the ON/OFF button on Isiris Monitor for at least 2 seconds. When Isiris Monitor is OFF restart by pressing ON/OFF button once more.
A recorded image is shown in the yellow file management tab.	Return to live image by pressing the blue Live image tab or the red button with a black X in the top right corner. Restart Isiris Monitor by pressing the ON/OFF button for at least 2 seconds. When Isiris Monitor is OFF restart by pressing ON/OFF button once more.
Isiris a is damaged.	Replace Isiris a with a new one.

Low picture quality.

Cause	Action
Light reflecting on Isiris Monitor screen.	Move Isiris Monitor to a position where no direct light influences the screen.
Dirty/damp screen.	Wipe the screen with a clean cloth.
Brightness and contrast settings not optimal.	Adjust the contrast and brightness using the designated menu on Isiris Monitor.
Blood, urine etc. on the lens (distal tip).	Gently rub the distal tip against the mucosa. If the lens cannot be cleaned this way remove Isiris a and wipe the lens with sterile gauze.

Absent or reduced irrigation capability.

Cause	Action
Leakage at Luer lock.	Reconnect properly.
Clamped irrigation line.	Unclamp.
Lack of pressure.	Elevate the bag.

Appendix 1: Electromagnetic compatibility

Like other electrical medical equipment the Isiris System requires special precautions to ensure electromagnetic compatibility with other electrical medical devices. To ensure electromagnetic compatibility (EMC) the Isiris System must be installed and operated according to the EMC information provided in this manual.

The Isiris System has been designed and tested to comply with IEC 60601-1-2 requirements for EMC with other devices.

Appendix 2. Warranty

There is no warranty on Isiris α.

STORE AWAY FROM LIGHT IN A COOL AND DRY PLACE.
DO NOT USE IF PACKAGE IS DAMAGED.
DO NOT RESTERILIZE.
DISCARD ISIRIS α AFTER USE.

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