

# SINGLE-USE LIKE YOU'VE NEVER SEEN IT BEFORE



Ambu® aScope™ 5 Broncho

**Ambu**



# SINGLE-USE BRONCHOSCOPY

Ambu® aScope™ 5 Broncho combines the familiar benefits of reusable scopes with the advantages of the single-use solution. It is compatible with the Ambu® aBox™ 2 displaying and processing unit (with full HD resolution), which can easily be connected to external monitors.

Ambu® aScope™ 5 Broncho

**BRINGS SINGLE-USE BENEFITS TO THE BRONCHOSCOPY SUITE**

## EXCELLENT MANOEUVRABILITY

### Access difficult-to-reach areas of the bronchial tree:


The combination of the bending angles, distal tip angulation range and flexibility of the insertion tube helps you access difficult-to-reach areas in the bronchial tree.

### 120° left/right rotation:

Designed to enable you to adjust the position of the distal end in order to ease navigation and insertion of endotherapy instruments.

### Broad endotherapy instrument compatibility:

Designed for use with the diagnostic and therapeutic endotherapy instruments commonly used in the bronchoscopy suite, including high-frequency tools.



## NEXT-LEVEL IMAGING

### High-resolution camera chip and full-HD displaying unit:

aScope 5 Broncho and aBox 2 work together to deliver superior imaging.

### Excellent image quality:

Advanced image processing, 120° field-of-view and 3-100 mm depth-of-field HD system with advanced adaptive auto control give you the bright and clear image needed for the procedure.

### Enhanced visibility:

Ambu's proprietary Advanced Red Contrast (ARC) enhancement algorithm is intended to improve the visibility of red colour tones in the image without altering the natural appearance.

### Easy to take, review and document images and videos:

At the touch of a button on the bronchoscope, you can zoom in and out, capture images and videos and review them post-procedure. Compatibility with PACS servers and EMR software makes documentation smooth and efficient.



## CONTACT US TO ARRANGE A DEMO

It's quick and simple for us to set up aScope 5 Broncho and aBox 2 at your place of work so you can get it in your hands and see for yourself.



Ambu® aScope™ 5 Broncho

# SINGLE-USE LIKE YOU'VE NEVER SEEN IT BEFORE

Sets a new standard in bronchoscopes

## 5 WAYS AMBU SINGLE-USE EFFICIENCY CAN BENEFIT WORKFLOW IN THE BRONCHOSCOPY SUITE

- **No repair hassle with single-use**

With the single-use aScope 5 Broncho solution, you get a new scope for every patient. This means there are no concerns about causing damage to aScope 5 Broncho with sharp or active endotherapy instruments.

- **Increase number of procedures without a sizable investment**

Freeing up costs and resources paves the way for your hospital to use the resources for other important tasks. With the aScope 5 Broncho solution, the only item needed to increase capacity is an additional displaying and processing unit.

- **Frees up costs and resources**

aScope 5 Broncho reduces the time and resources spent on cleaning. In fact, it completely eliminates the need for reprocessing. Along with that, the cost and effort of maintenance and repairs, service contracts, and resources spent on transportation, administration, training and cleaning supplies are removed.

- **Easier to plan and manage your schedule**

With the sterile, single-use aScope 5 Broncho, you have a variety of sizes available in storage whenever you need them. This can save you from the inefficiency and bother of cancellations and rescheduling.

- **An extra measure of safety**

The aScope 5 Broncho solution eliminates the risk of patient-to-patient contamination. This could be especially relevant in reducing the risk of infection transmission among immunosuppressed pulmonary patient populations.

**Ambu** FOREVER FORWARD

Ambu A/S  
Baltorpbakken 13  
DK - 2750 Ballerup  
Denmark  
Tel.: +45 7225 2000  
ambu.com