

SINGLE-USE LIKE YOU'VE NEVER SEEN IT BEFORE

A NEW ERA FOR SINGLE-USE BRONCHOSCOPY



Ambu® aScope™ 5 Broncho

Ambu



ONE SOLUTION FOR A BROAD RANGE OF PROCEDURES

With aScope 5 Broncho, you're not limited by the number of available endoscopes. You have a variety of sterile bronchoscopes in different sizes available whenever you need them for diagnostic, therapeutic and peripheral bronchoscopy procedures.

Diagnostic and therapeutic procedures

The aScope 5 Broncho HD 5.0/2.2 and 5.6/2.8 endoscopes offer excellent manoeuvrability, even when working with endotherapy instruments, including high-frequency tools. This, combined with a bright and clear image, makes them suitable for a wide range of diagnostic and therapeutic procedures.



BAL procedures

The single-use Ambu® aScope™ 5 Broncho Sampler Set is a unique integrated sampling solution that protects the sample, simplifies workflow and improves safety.

It includes a 60 ml sample container to address the need for larger samples in BAL procedures.



The Ambu® aScope™ 5 Broncho Sampler Set

Protects the sample from start to finish

- Closed-loop system reduces the risk of sample contamination and loss common with Lukens traps
- Minimises exposure points
- Protects the HCP from exposure to airway secretions

Simple to operate without assistance

- Seamlessly integrated with the aScope 5 Broncho endoscope
- Allows you to switch between suction and sampling at the turn of a knob, without the need to switch tubing
- High degree of control when navigating the lungs

Sterile and ready to go

- Single-use solution is sterile straight from the package
- Always available when you need it
- Clicks right onto the aScope 5 Broncho endoscope



Peripheral bronchoscopy procedures

aScope 5 Broncho 4.2/2.2 offers a unique opportunity to take advantage of a regular-sized working channel in a thin bronchoscope. The thin size and high bending angles of 210°/210° enable excellent manoeuvrability.

What's more, it is compatible with commonly used endotherapy instruments and active tools*. These capabilities make aScope 5 Broncho 4.2/2.2 ideal for peripheral bronchoscopy procedures.

Smaller anatomies

With its thin outer diameter, aScope 5 Broncho 2.7/1.2 can be used in smaller anatomies without compromising on manoeuvrability, image quality or access to the working channel. In fact, it offers a 120° field of view and 3-100mm depth of field. The ultrathin aScope 5 Broncho 2.7/1.2 is suitable for ET tubes as small as 3.5mm.



STERILITY AND EFFICIENCY OF THE SINGLE-USE CONCEPT

- **Sterile:** No risk of patient-to-patient contamination
- **Always available:** No reprocessing or repair bottlenecks
- **Consistent performance:** No wear and tear affecting the bending performance or imaging



THE LEVEL OF QUALITY NEEDED IN THE BRONCHOSCOPY SUITE

- **High performance:** Excellent manoeuvrability and high-quality imaging
- **Greater flexibility:** A variety of readily available bronchoscopes in 4 different sizes
- **A unique sampling solution:** The aScope 5 Broncho Sampler Set simplifies BAL procedures (see other side)

*aScope 5 Broncho 2.7/1.2 and 4.2/2.2 are not compatible with high-frequency tools

Ambu® aScope™ 5 Broncho

SINGLE-USE LIKE YOU'VE NEVER SEEN IT BEFORE

A new era for single-use bronchoscopy

aScope™ 5 Broncho flexible single-use bronchoscope

Item code	Product	Quantity
622001000	Ambu® aScope™ 5 Broncho HD 5.6/2.8	5
621001000	Ambu® aScope™ 5 Broncho HD 5.0/2.2	5
620001000	Ambu® aScope™ 5 Broncho 4.2/2.2	5
624001000	Ambu® aScope™ 5 Broncho 2.7/1.2	5

aScope™ 5 Broncho Sampler Set flexible single-use bronchoscope and specimen sampling system

Item code	Product	Quantity
622002000	Ambu® aScope™ 5 Broncho HD 5.6/2.8 Sampler Set	5
621002000	Ambu® aScope™ 5 Broncho HD 5.0/2.2 Sampler Set	5
620002000	Ambu® aScope™ 5 Broncho 4.2/2.2 Sampler Set	5

Single-use sample container

Item code	Product	Quantity
480 500 000	Ambu® aScope™ BronchoSampler 60 SC	10

Displaying and processing units

Item code	Product	Quantity
505001000*	Ambu® aBox™ 2	1
405011000*	Ambu® aView™ 2 Advance	1

aCart™

Item code	Product	Quantity
9.713.800	Ambu® aCart™ Compact KU	1
3.000.991	Ambu® aCart™ Compact Vario Mount HA	1

*Item code differs by region. Please refer to local material or contact Ambu to obtain the full item code for placing orders.

Ambu FOREVER FORWARD

Ambu A/S
Baltorpbakken 13
DK - 2750 Ballerup
Denmark
Tel.: +45 7225 2000
ambu.com